Temperature

Controlled Shaking

Leading Safety Standards

Superior Ease of Use

Reduced Cost of Ownership

Operating Manual

Page 2

Platform Shakers Hei-MIX Unimax 2010/Promax 2020/
Polymax 2040

Operating manual must be read before initial start-up.
Please follow the safety instructions provided.
Please keep for future reference.

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About this Document	4
Version and Variants	4
About this manual	4
Terms and signs used.....	5
Safety Instructions	6
General safety instructions.....	6
EC Declaration of Conformity.....	6
Intended use.....	6
Installation / Electrical safety	7
Personnel qualification.....	7
Operating company's obligations.....	7
Installation site	7
Changes to the unit	7
Safety for the personnel	8
Safety during use	8
Disposal	8
Device Description	9
Device overview	9
Hei-MIX Polymax 2040 overall view.....	9
Hei-MIX Unimax 2010 overall view	9
Hei-MIX Promax 2020 overall view	10
Hei-MIX Polymax 2040, Unimax 2010, Promax 2020 control panel	10
Start-Up	11
Setting up device	11
Connect/disconnect power cord	11
Switch unit on and off	12
Operation	13
Mounting vessels.....	13
Secure vessels on shaker plate	13
Shaking mode.....	13
Continuous mode.....	14
Timer mode	14

Cleaning and Maintenance	15
Cleaning.....	15
Maintenance	15
Troubleshooting	15
Possible errors and how to resolve	15
Assembly	16
Electrical connections.....	16
Mounting attachments.....	17
Base frame and tension rollers	17
Clamps for Erlenmeyer flasks	17
Clamps for separatory funnel.....	18
Multi-tier attachment	19
Disassembly and Storage	21
Dismantling, transportation and storage	21
Dismantling	21
Transportation and storage	21
Accessories and Spare Parts	22
Scope of delivery.....	22
Accessories.....	22
Attachments	23
Technical data	23
Performance range / admissible load	24
Service	25
Contact / Technical Service	25
Warranty	26
Confirmation of condition of unit.....	27

➤ Version and Variants

Version

This manual describes the function, operation and maintenance of the Platform Shakers Hei-MIX Unimax 2010/Promax 2020/Polymax 2040.

Version	Alteration date
1.3	07.2015

Variants

The units are available in different configurations. Certain features and functions apply to specific configurations only. The varieties are documented in this manual.

➤ About this manual

This operating manual is a component of the device described here.

- Please read this manual carefully and obey all safety and warning notices.
- Ensure that every operator reads this manual.
- Ensure that this manual is accessible for every operator.
- Pass on the operating manual to the subsequent owner.

For the current version of this operating manual in pdf format, please go to:

<http://www.heidolph-instruments.com/Support/Operation Manuals/Shakers & Mixers/Platform Shakers/Unimax 2010/Promax 2020/Polymax 2040>

In addition please observe the regional regulations.

➤ Terms and signs used

In this manual you will find standardized terms and signs to warn you of possible dangers or give important hints. Please take special notice to these signs to avoid personal injury or damage to the unit.

The following terms and type of signs (basic drawing elements) are used:

Symbol	Additional term / Description
Warning signs 	<p>The yellow triangle indicates hazardous situations. It is used in combination with the following terms:</p> <p>DANGER: Indicates a hazardous situation which, if not avoided, will result in irreversible injury or death.</p> <p>WARNING: Indicates a hazardous situation which, if not avoided, could result in serious injury or material damage.</p> <p>CAUTION: Indicates a hazardous situation which, if not avoided, may result in a minor injury or material damage.</p>
Prohibitory signs 	<p>Not to do:</p> <p>The red circle indicates a situation that should be avoided under all circumstances to avoid injury or damage.</p>
Mandatory signs 	<p>To do:</p> <p>The blue circle indicates important information. Please obey to avoid property damage.</p>

Other signs used

Symbol	Description
➔	Handling instruction, action required
✓	Result of action
▪	List of information
▪	
▪	
a. b. c.	List of variants

➤ General safety instructions

EC Declaration of Conformity

This device complies with the following EC-Directives*:

- 2006/42/EC Machinery Directive
- 2006/95/EC Low-Voltage Directive
- 2004/108/EC Electromagnetic Compatibility Directive

*Also see attachment „EC Declaration of Conformity“.

The device has been constructed according to state-of-the-art technology and recognized safety regulations. However, risks may still arise during installation, operation and maintenance.

→ Please ensure the operating manual is available at all times.

The device may only be used under the following circumstances:

- Only operate the device, if it is in full working order.
- Ensure all operators of the device possess the necessary safety and risk awareness.
- Operate the device according to the instructions stipulated in this manual only.
- If there is something you do not understand, or certain information is missing, ask your manager or contact the manufacturer.
- Do not do anything on the device without authorization.
- Only use the device according to its intended use.

Intended use

The device is intended for use by trained and authorized personnel only.

The device is suitable for the following use:

- mixing
- shaking
- separating
- blending of solvents

The operation for the device is suitable in one of the following locations:

- in any research application

Improper use

Any use which deviates from the device's intended use is considered to be improper. The manufacturer does not accept liability for any damages resulting from non-permitted uses. The risk is carried by the operator alone.

Installation / Electrical safety

- The device may only be connected when the mains voltage corresponds to the information on the type plate of the unit.
- The mains connection must be accessible at all times.
- Repairs may only be performed by a qualified electrician.
- Never operate the unit with a damaged power cord.
- Always turn the unit OFF and disconnect mains power before performing any maintenance or service.

Personnel qualification

- The device may only be operated by qualified persons.
- The device may only be operated by individuals who have been instructed in its proper use by qualified persons.
- The device may only be operated and maintained by persons who are of legal age.
- Other personnel may only work with the unit under continuous supervision of an experienced qualified operator.
- This manual must be read and understood by all persons working with the device.
- The personnel must have received special safety instructions in order to guarantee responsible and safe work procedures.

➤ Operating company's obligations

Installation site

- The device must be positioned in a suitable location.
- The device must be installed sufficiently stable.
- All screw connections must be securely tightened.
- Installation and operation of the device is only permitted in facilities which are fitted with the appropriate laboratory equipment (e.g. with air extraction units).
- We recommend operation in fume hood enclosures.
- The device may only be operated in enclosed spaces and under the following environmental conditions:

Ambient temperature

- 5 - 40 °C
 - 80 % relative humidity, no condensation at 5 - 31 °C, linear decreasing to 50 % at 40 °C
-

Changes to the unit

- No unauthorized changes may be made to the unit.
- No parts may be used which have not been approved by the manufacturer.
- Unauthorized changes result in the EC Declaration of Conformity losing its validity, and the unit may no longer be operated.
- The manufacturer is not liable for any damage, danger or injuries that result from unauthorized changes or from operating the unit other than described in this manual.

Safety for the personnel

- Ensure that only qualified personnel operates the device.
- Observe the following regulations:
 - Laboratory guidelines
 - Accident prevention regulation
 - Ordinance on Hazardous Substances
 - Other generally accepted rules of safety engineering and occupational health
 - Local regulations

➤ Safety during use

- Wear the appropriate protective clothing when working on the device (clothing, protective glasses and if necessary, safety gloves).
- Do not use the device in potentially explosive areas. The device is not protected against explosion. There is no explosion or ATEX protection available.
- Do not operate or assemble devices in the vicinity which are emission or radiation sources (electromagnetic waves) for the frequency range (3×10^{11} Hz to 3×10^{15} Hz).
- Avoid putting pressure on the display when you do not operate the device.
- Eliminate errors immediately.
- Do not use abrasive material to clean the glass surfaces. Only wipe with damp cloths.
- Always switch the device off after use.

➤ Disposal

- Check the device components for hazardous substances and solvents.
- Clean all components before disposal.
- Dispose of the device according to the valid national legal regulations.
- Dispose of the packaging material in accordance to the valid national legal regulations.

› Device overview

Hei-MIX Polymax 2040 overall view

Hei-MIX Unimax 2010 overall view

Start-Up

Hei-MIX Promax 2020 overall view

Hei-MIX Polymax 2040, Unimax 2010, Promax 2020 control panel

➤ Setting up device

Warning:
Risk of slipping device!

As a result of vibrations generated during operation the device might slide off the table top.

The surface must be smooth and clean and should be properly leveled.

- Locate the shaker on a stable, horizontal surface.
- Clean surface and feet with a damp cloth and ethanol regularly.

➤ Connect/disconnect power cord

Connect power cord

The power cord comes with a three-pole plug and a recess on the bottom side.

Unit plug

Power cord

The socket is located on the back side of the unit.

- Connect plug to the socket with the recess facing down.
- Push plug up against socket until securely attached.
- Connect plug to power socket.

Unit socket

Disconnect power cord

- Disconnect plug from power socket first.
- Then disconnect plug from socket on unit.

➤ Switch unit on and off

The ON/OFF switch is situated at the front panel.


Switch unit on

- Press the On/Off switch at the lower side.
- ✓ The green control light switches on.
- ✓ The unit is switched on.


Switch unit off

- Press the On/Off switch at the upper side.
- ✓ The green control light switches off.
- ✓ The unit is switched off.

Toggle switch ON/OFF

Unit switched on

Unit switched off

➤ Mounting vessels

Secure vessels on shaker plate

Adaptors for Erlenmeyer flasks, separating funnels and an adaptor with clamping rollers for universal use are available.

- To secure vessels on shaker plate, use optional adaptors (see chapter „Spare parts and accessories“ and chapter „Assembly“).
- At low shaking speed, low-profile vessels like petri dishes may be placed directly on the non-skid rubber plate.

A single flask should be arranged in the middle of the shaker plate.
Several flasks should be distributed equally on the plate.

➤ Shaking mode

Caution: **Risk of breakage and spillage**

If the vessels are not attached securely they may fall out and spill the sample.

Prior to switching the unit on ensure that all vessels are securely fastened.

If working with separatory funnels make sure that the outlet pipes of your separatory funnels are clear of other items on your lab bench.

Warning: **Risk of poisoning**

Open vessels and too high vibration frequencies may result in samples splashing.

Close all vessels tight especially if sample is dangerous or toxic.

Adjust vibration frequency step by step until you have reached the required rpm settings..

Wear safety glasses and especially with dangerous and toxic samples appropriate protective clothing.

At high loads (load-bearing capacity of shaker plate) and high frequency of shaking, always make reference to load graph as applicable for your item (see chapter „Technical data“).

When using shakers inside conditioning cabinets, make reference to ambient conditions as stipulated in chapter „Technical data“.

Operation

You can choose from two shaking options:

- Continuous mode: The unit starts and stops shaking as soon as you switch it on or off.
- Timer mode: The time for shaking can be set between 0 and 120 minutes. The unit stops shaking when the time has elapsed.

Continuous mode

- ✓ Select mode switch is positioned to continuous mode „0”.
- Turn item on with ON/OFF switch.
- Set speed with analog knob „Speed”.
- ✓ Actual speed appears in the 3-digit digital display.
- Turn speed off.
- Turn item off with ON/OFF switch.

Analog knob „Speed” Analog knob „Time”

Select mode switch:
continuous / timer

Timer mode

- Select timer mode „**Timer**” with select mode switch.
- Set time with analog knob „Timer”.
- Set speed with analog knob „Speed”.
- Turn item on with ON/OFF switch.
- ✓ When time has elapsed an acoustic signal will be audible and shaking action will stop.
- Turn speed off.
- Turn item off with ON/OFF switch.

A power cut will stop shaking mode. The timer will continue running even in case of power supply failure.

> Cleaning

All surfaces of the unit allow for cleaning with a damp cloth and if required a mild soap lotion.

Attention: **Damage to the surfaces**

All surfaces can be damaged by improper cleaning.

In any case do not use:

- Chlorine bleach or any cleaning agent including chlorine
- Ammonia
- Abrasive cleaning agents such as cleaning rags, scrubbing agents or any other agents which include metal components

As a result of vibrations generated during operation the device might slide off the table top.

To ensure a stable position clean surface and feet with a damp cloth and ethanol regularly.

> Maintenance

The unit is maintenance-free. Any necessary repair must be performed by an authorized Heidolph distributor.

Please contact Heidolph Instruments or your local Heidolph Instruments distributor.

> Possible errors and how to resolve

Malfunctions	Possible reason	Troubleshooting
On-/Off switch does not enlighten	No power	<ul style="list-style-type: none"> ▪ Check size and shape of plug and compatibility to your electrical socket ▪ Check circuit breakers
	Light diode defective	Contact your local Heidolph Instruments distributor
Shaking movement stops	Pre-set timer elapsed	Restart timer, if necessary
	Thermal motor circuit breaker triggered by motor overload.	Wait 20 minutes, decrease overall load
	<ul style="list-style-type: none"> ▪ Motor defective ▪ Control electronics defective 	Contact your local Heidolph Instruments distributor

If you experience a malfunction which can not be resolved, please contact your authorized Heidolph Instruments distributor immediately.

➤ Electrical connections

- This unit must **only** be connected **by a professional electrician**.
- **Any repair** on this unit must **only** be carried out **by a qualified specialist**.

Warning: Electric shock!

This unit must only be connected to a grounded electrical socket.

Repairs are allowed by qualified professionals only.
Any improper repair can result in a dangerous situation. Contact your local Heidolph Instruments distributor for any repair you may have.

The unit and voltage must match. The type label on the back side of the unit provides all voltage specifics.

When shipped, the unit is grounded. In case the original plug is changed, the new plug must have a protective conductor!

Rating plate

Fuse

Unit socket

Color code for electrical connections:

Green/Yellow	PE: Protective conductor (E arth)
Blue	N: Neutral conductor
Brown	P: Phase
Black	P: Phase
Grey	P: Phase

The unit is connected with the power cord supplied to the electrical socket. It plugs into the housing on the rear side of the unit.

For countries utilizing a plug other than the standard one supplied:

- The plug supplied may only be changed by a professional electrician.
- If utilizing an adaptor ensure that it is approved by the local regulations.

➤ Mounting attachments

Base frame and tension rollers


Base frame and tension rollers allow to safely attach whatever sort of vessels and racks on the shaker platform.

Components frame tension roller:

- 4 bolts
- 1 box-end spanner
- 2 holding bars

Necessary on top:

- Min. 2 pcs tension rollers


Bolts

Holding bars

Tension rollers

- ➔ Thread 4 ea. bolts in the respective holes of platform and pull tight with box-end spanner from hardware bag.
- ➔ Slip holding bar over bolts and secure with wing bolts.
- ➔ Slip tension rollers over holding bars, push against the vessel to be clamped, and secure with wing bolts.
- ➔ Attaching one more row of vessels requires installation of extra tension rollers.
- ➔ Lock tension rollers safely in place and pull all screws of the universal mount tight as specified.


Clamps for Erlenmeyer flasks

The perforated platform with appropriate clamps allows to safely hold Erlenmeyer flasks of even different sizes on the shaker plate.

Distribution of holes in the rack permits whatever combination of clamps.


- ✓ At the bottom of the perforated platform a box-end spanner is attached to mount clamps.


Assembly

- 25 to 100 ml clamps feature **one** threaded bolt.
 - 250 - 1000 ml clamps feature **four** threaded bolts.
- Remove nuts and washers from clamp bolts.

Nuts Washers Clamp bolts


- Plug clamp in rack hole.
- Secure with washer.
- Tighten nuts with box-end spanner.
- Ensure all screws have been pulled tight.

Nuts Washers Clamp bolts

The only restriction is to arrange clamps in a manner avoiding direct contact between each other; moreover, flasks should not protrude over the edge of the shaker plate.

- Locate perforated platform with clamps installed on shaker plate.
- Make sure, crimped edge is plugged over shaker plate, thus positively locking both parts with each other.

Lightly tilt Erlenmeyer flasks to fit them under the clamp.

Clamps for separatory funnel

Perforated platform and separatory funnel clamp safely hold size 250 ml to 2000 ml separating funnels on the shaker plate.

The perforated platform accepts the following number of separatory funnel clamps:

- 4 ea. 250 ml
- 3 ea. 500 ml
- 3 ea. 1000 ml
- 2 ea. 2000 ml

A mix of various sizes may be used either.

Assembly


Caution Risk of breakage

If clamps for separatory funnels are mounted too close, separatory funnels could break when touching each other during operation.

Maintain sufficient space between separating funnels to avoid a collision situation.

- Remove nuts and washers from side bars B and C.
- Plug side bars B and C into the perforated platform holes, so that dimension „A“ is about 25 mm less than length „L“ of the separatory funnel.
- ✓ The distance should guarantee proper clamping of the separatory funnel and its closing plug by means of compression spring.
- ✓ Dimension „A“ once set will be good for one and the same size of separatory funnel only.
- ✓ Whenever using a different model of separating funnel, you need to check dimension „A“ and readjust accordingly.
- Install washer and nuts on bolts of side bars B and C; pull tight with box-end spanner from your hardware bag.
- Locate perforated platform with funnel clamps on shaker plate.
- Make sure, crimped edge is plugged over the shaker plate, so positively locking both parts with each other

Install separatory funnel

- Install separating funnel with its outlet pipe first through compression spring.
- Press item against compression spring until closing plug locks with retainer ring.

Multi-tier attachment

The multi-tier attachment allows installation of a second or third rack plate for Erlenmeyer flasks on the shaker plate.

- 25 ml Erlenmeyer flasks: allows a total of **3** perforated platforms
- 50 ml bis 250 ml Erlenmeyer flasks: allows a total of **2** perforated platforms



Assembly

Components multi-tier attachment:

- 4 columns
- 1 box-end spanner
- 8 clamping sleeves

To adjust height of clamping sleeves correctly the columns feature four marks A, B, C and D.

- 25 ml Erlenmeyer flask, 3 platforms, marks A and D
 - 50 ml Erlenmeyer flask, 2 platforms, marks C
 - 100 ml Erlenmeyer flask, 2 platforms, marks B
 - 250 ml Erlenmeyer flask, 2 platforms, marks A
- Install 4 ea. columns and pull tight with box-end spanner from your hardware bag.
- Locate first platform over columns.
- Slip 4 ea. clamping sleeves with wing bolt situated below over columns.
- Align top rim of clamping sleeves with adequate marks and secure with wing bolts.
- ✓ Wing bolts must be adjusted vertically to the outer edge of the platform.
- Ensure all bolts and screws have been pulled tight.
- Install second platform over columns and seat on clamping sleeves.
- Install 3rd platform (if required) by repeating above steps.

➤ Dismantling, transportation and storage

Dismantling

Warning:
Risk of injury!

While the unit is plugged in you might be at harm if the device is accidentally switched to ON mode.

Prior to dismantling always disconnect the plug from the mains.

- ✓ The unit is switched off.
- Disconnect plug from power socket.
- Disconnect plug from socket on the unit.
- Remove all vessels.
- Disassemble attachments.

Transportation and storage

- Store and transport the unit and its components only if they were emptied and cleaned in the original packing material.
- Alternatively store and transport the unit in a container which prevents damages during transportation.
- Seal the packing carefully against unauthorized or mistaken opening.
- Store the unit in a dry place.

Improper transportation may result in damages to the system and the mechanical components.

Avoid any kind of shocks during transportation.

➤ Scope of delivery

Component	Variant	Quantity	Product number
Hei-MIX	Unimax 2010	1	542-10020-00
	Promax 2020	1	542-20020-00
	Polymax 2040 tilt angle 5°	1	542-40005-00
	Polymax 2040 tilt angle 10°	1	542-40010-00
Operating Manual English / German		1	01-005-005-68
Warranty registration / Confirmation of condition		1	01-006-002-78
EC Declaration of Conformity		1	01-001-025-06
Power cord		1	

➤ Accessories

Component	Quantity	Product number
Frame tension roller (without tension roller)	1	549-50000-00
Tension roller	1	549-58000-00
Universal perforated platform 2000	1	549-59000-00
Multi-tier attachment	1	549-62000-00
Clamps for Erlenmeyer flasks, sizes from 25 ml - 2000 ml		See general catalog

Additional accessories are shown in the general catalog or on our website at www.heidolph-instruments.com.

➤ Technical data


Supply power (W)	1 x 230 V (50/60 Hz)
Operational voltage	100 - 240 V (50/60 Hz)
Power consumption (W)	115
Protection class (DIN EN 60529)	IP 20
Drive	EC-motor
Overheat protection	self-resetting
Operating mode	timer or continuous
Speed setting	digital
Platform size (B x T)	390 x 340
Weight (kg)	16
Load capacity (kg)	10
Permissible ambient conditions	<ul style="list-style-type: none"> ▪ 5 - 40 °C ▪ 80 % relative humidity, no condensation at 5 - 31 °C, linear decreasing to 50 % at 40 °C

Model Hei-MIX	Unimax 2010	Promax 2020	Polymax 2040
Motion	orbital	reciprocating	wave
Dimensions (w x h x d) (mm)	426 x 135 x 435		426 x 208 x 435
Speed range (rpm)	20 - 400		2.5 - 50
Orbit / Stroke (mm)	20		-
Tilt angle (°)	-		5 / 10


Performance range / admissible load


Hei-MIX Unimax 2010

Hei-MIX Promax 2020

Hei-MIX Polymax 2040

➤ Contact / Technical Service

Questions / Repair work

If any aspect of installation, operation or maintenance remains unanswered in the present manual, please get in touch with the following address.

For repairs please call Heidolph Instruments or your local authorized Heidolph Instruments distributor.

Warning: Danger of poisoning!

Contaminated units can lead to severe injury or death of our employees!

When shipping items for repair that may have been contaminated by hazardous substances, please:

- advise exact substance
 - take adequate protective action towards our parts receiving and service personnel
 - mark the pack in accordance with Ordinance on Hazardous Substances
-

Included at the end of this manual you will find a „Confirmation of condition of unit“.

→ Kindly copy and fill in this form and submit it prior to shipping the unit for repair.

Contact details:

Heidolph Instruments Germany

Heidolph Instruments GmbH & Co. KG
Technical Service
Walpersdorfer Str. 12
D-91126 Schwabach / Germany

Tel.: +49 - 9122 - 9920-74
Fax: +49 - 9122 - 9920-84

E-Mail: service@heidolph.de
www.heidolph.de

Heidolph Instruments North America

Phone: 1-866-650-9604
E-mail: service@heidolph.com
www.heidolphNA.com

Heidolph Instruments United Kingdom

Phone: 01799 - 5133-20
E-mail: service@radleys.co.uk
www.heidolph-instruments.co.uk

All other countries

You will find contact details of your local Heidolph distributor at
www.heidolph.com

➤ Warranty

Heidolph Instruments provides a 3-year warranty for the products described here (excluding glass and wearing parts) if you register using the warranty card enclosed or online (www.heidolph.com). The warranty is valid from the point of registration. The serial number is also valid without registering.

The warranty covers part and manufacturing defects.

In the case of a part or manufacturing defect, the device shall be repaired or replaced free of charge under the terms of the warranty.

Heidolph Instruments shall not assume liability for any damages incurred as a result of improper handling or transport.

Warranty claim?

- ➔ Please inform Heidolph Instruments or your local Heidolph distributor should you wish to make a warranty claim.

➤ Confirmation of condition of unit

→ In the case of repair, copy and complete the Confirmation of condition of unit and send it to Heidolph Instruments.

1. Details about the unit

Product number

Serial number

Reason for repair

2. Has the device been cleaned, decontaminated/sterilized?

Yes

No

3. Is the unit in a condition which does not represent any health threats for the staff of our service department?

Yes

No

If not, which substances has the unit
come into contact with?

4. Legally binding declaration

The customer is aware of being legally liable to Heidolph
Instruments for any damages arising from incomplete and
incorrect information.

Date

Signature

Company stamp

Please note

The shipper is responsible for the return of the goods in well
packed condition, suitable for the mode of transport.

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This is a translation of the original operating manual. We reserve the right to make any technical modifications without prior notification. This documentation is not subject to revision; the latest version of this documentation can be found on our homepage. The copyright for all texts and images is held by Heidolph Instruments GmbH & Co. KG.
© Heidolph Instruments GmbH & Co. KG